



Rotavirus & Adenovirus Antigen Rapid Test Kit (Colloidal Gold Method)

The Rotavirus & Adenovirus Antigen Rapid Test Kit is a lateral flow immunoassay assay designed for qualitative detection of rotavirus and adenovirus antigens in human fecal specimen. This method is easy to operate, and the result can be obtained by observing whether there is a red and purple strip in 15-20 minutes. Simple operation steps. It is intended to be used as a screening test and provides a preliminary test result to aid in the diagnosis of infection with rotavirus and adenovirus.

Test Procedure

- **Specimen:** Whole blood, serum or plasma.
- **Specimen collection:** Specimens of all types are collected by conventional methods.
- **Test method:**



Specifications

Components	Test card/ Sample extraction tubes / Sample Droppers
Pack Size	1Test/Kit;25Tests/Kit; 50 Tests/Kit
Storage Conditions	2°C~30°C; Valid for 24 months
Test samples	Fecal
Reaction conditions	15~30°C 15~20 minutes
Special instrument	No need

Product Informations

Product Name	Cat.No	Package	Sample Diluent	Sample Droppers	Sample Type
Rotavirus & Adenovirus Antigen Rapid Test Kit (Colloidal Gold Method)	BSK17T1S	1 Test/Kit	1 mL/Tube	1	Fecal
	BSK17S1S	25 Tests/Kit	25 × 1 mL/Tube	25	
	BSK17M1S	50 Tests/Kit	50 × 1 mL/Tube	50	